



DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set

Instructions for use

Version: 1.0

Research Use
Only

Qingdao MGI Tech Co., Ltd.

About the instructions for use

This instructions for use is applicable to DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (hereinafter referred to as preparation set). The instructions for use version is 1.0, and the preparation set version is V1.0.

This instructions for use and the information contained herein are proprietary to Qingdao MGI Tech Co., Ltd. (hereinafter referred to as MGI), and are intended solely for the contractual use of its customers in the use of the products described herein and for no other purpose. This instructions for use and its contents shall not be reprinted, reproduced, modified, distributed, or disclosed to others, in whole or in part, without prior written consent from MGI.

MGI makes no commitment to this instructions for use, including (but not limited to) any special commercial purpose and any reasonable implied warranties. MGI has taken measures to ensure the correctness of this instructions for use. However, MGI is not responsible for any missing parts of the instructions for use, and reserves the right to revise the instructions for use and modify the preparation set, so as to improve the reliability, performance or design.

Figures in this manual are for illustrative purpose only. The content may be slightly different from the device. For the most up-to-date details, refer to the preparation set purchased.

Qubit™ are trademarks of Thermo Fisher Scientific, or its subsidiaries. DNBSEQ™ is a trademark of MGI or its subsidiaries in China and/or other countries. Other names and trademarks mentioned in this instructions for use may be the property of their respective companies and subsidiaries.

©2023 Qingdao MGI Tech Co., Ltd. All rights reserved.

Revision history

	Date	Version
Initial release	July 3, 2023	1.0



Tips

- Please download the latest instructions for use, and use it with the corresponding version of the preparation set.
- Download the instructions for use through search according to the catalog number or product name from the website:
www.mgi-tech.com/download/files.

Contents

Chapter 1 Product overview	1
1.1 About the product	1
1.2 Intended use	1
1.3 Compatible sequencing platform	1
1.4 Component information	2
1.5 Storage and transportation condition	4
1.6 Self-provided materials	5
1.7 Precautions	7

Chapter 2 Sample requirements and processing	8
2.1 Precautions before experiment	8
2.2 Pre-experiment preparation	8
2.2.1 Sample requirements	8
2.2.2 Experiment requirements	9
2.2.3 Preparing reagents	9
2.3 Preparing nuclei suspension	10

Chapter 3 Nuclei transposition	11
3.1 Pre-experiment preparation	11
3.2 Nuclei transposition	11

Chapter 4 Droplet formation	12
4.1 Pre-experiment preparation	13
4.1.1 Preparing reagents and consumables	13
4.1.2 Preparing sample phase suspension	13
4.1.3 Preparing beads phase suspension	14
4.2 Performing droplet formation	15

Chapter 5 Performing PCR in droplets	18
---	-----------

5.1 Pre-experiment preparation	19
5.2 Recovering droplets	19
Chapter 6 Demulsification and beads selection	20
6.1 Pre-experiment preparation	21
6.2 Performing demulsification	21
6.3 Beads selection	22
Chapter 7 Library preparation	23
7.1 Pre-experiment preparation	23
7.2 Performing library preparation	24
7.3 Beads selection	25
Chapter 8 Circularization	27
8.1 Pre-experiment preparation	27
8.2 Denaturation	27
8.3 Single-stranded circularization	28
8.4 Enzyme digestion	29
8.5 Beads purification	29
Chapter 9 Sequencing	31
9.1 Sequencing requirements of DNBSEQ-G400RS	31
9.1.1 Pre-experiment preparation	31
9.1.2 Making DNB	32
9.1.3 Pooling libraries	32
9.1.4 Sequencing parameters	32
9.2 Sequencing requirements of DNBSEQ-T7RS	34
9.2.1 Pre-experiment preparation	34
9.2.2 Making DNB	35
9.2.3 Pooling libraries	35
9.2.4 Sequencing parameters	35

Appendix 1 About the DNA Clean Beads and purification	36
<hr/>	
Appendix 2 Using scATAC Barcode Primer (1 to 16)	37
<hr/>	
Appendix 3 96 sequences information and concentration of scATAC Barcode Primer (1 to 96)	40
<hr/>	
Appendix 4 Manufacturer information	47

---This page is intentionally left blank.---

Chapter 1 Product overview

This chapter describes basic information of the product, including the intended use, compatible sequencing platforms, component information, and so on.

1.1 About the product

Based on the proprietary DNBelab C series single-cell library preparation and DNBSEQ™ sequencing technologies, combined with self-developed single-cell analysis software, the DNBelab C Series Cell Omics Solution can realize a portable, instant, and one-stop single-cell omics research workflow.

Based on droplet microfluidics, the DNBelab C Series High-throughput single-cell ATAC library preparation set, to quickly process the single-cell or single-cell nucleus suspension and prepare the dedicated libraries applicable to DNBSEQ sequencing platforms of MGI by catching beads efficiently. All reagents, chips, and consumables included in the preparation set have undergone strict quality control and function verification, to ensure single-cell ATAC library preparation stability and reproducibility.

1.2 Intended use

This preparation set is applicable to preparation of high-throughput single-cell ATAC libraries from human or mammal animal samples. Before use, it is necessary to prepare the samples into single-cell nuclei suspensions.



WARNING This preparation set is for scientific research use only, and cannot be used for clinical diagnosis.

1.3 Compatible sequencing platform

Sequencing platform	DNBSEQ-G400RS
	DNBSEQ-T7RS
Sequencing recipe	115 (Read 1) + 69 (Read 2) + 10

1.4 Component information

The preparation set includes the following 3 reagent kits and 1 chip kit.

Table 1 DNBelab C Series High-throughput Single-cell ATAC Library Preparation set (16 RXN) (Cat. No.: 940-000793-00)

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 1 Droplet Formation Kit) (Cat. No.: 940-000792-00)	scATAC Cell Beads	Black	1.6 mL/tube × 1
	scATAC Wash Buffer	Brown	1.6 mL/tube × 2
	scATAC Lysis buffer	Black	768 µL/tube × 1
	scATAC Breakage Reagent	Brown	1.6 mL/tube × 1
	P100 Oil	Natural	11.2 mL/bottle × 1
	Mineral Oil	Natural	6.4 mL/bottle × 1
	DNA Clean Beads	Natural	5.933 mL/bottle × 1
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 2 Droplet Formation Kit) (Cat. No.: 940-000794-00)	Tn5 Transposase	Green	64 µL/tube × 1
	5 × TAG Buffer	Green	80 µL/tube × 1
	scATAC Bead Buffer	Black	1.114 mL/tube × 1
	PCR Primer	Black	39 µL/tube × 1
	PCR Enzyme I	White	128 µL/tube × 1
	scATAC Cell Buffer	White	748 µL/tube × 1
	scATAC DIR Reagent	White	32 µL/tube × 1
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 3 Library Preparation Kit) (Cat. No.: 940-000910-00)	PCR Enzyme II	Blue	800 µL/tube × 1
	scATAC Splint Buffer	Yellow	80 µL/tube × 1
	scATAC Cyclization Buffer	Yellow	157 µL/tube × 1
	scATAC DNA Ligase	Yellow	8 µL/tube × 1
	scATAC Exo Enzyme	Orange	42 µL/tube × 1
	scATAC Exo Buffer	Orange	23 µL/tube × 1
	scATAC Stop Buffer	Orange	48 µL/tube × 1
	scATAC Barcode Primer 1	Red	16 µL/tube × 1
	scATAC Barcode Primer 2	Red	16 µL/tube × 1
	scATAC Barcode Primer 3	Red	16 µL/tube × 1

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 3 Library Preparation Kit) (Cat. No.: 940-000910-00)	scATAC Barcode Primer 4	Red	16 µL/tube × 1
	scATAC Barcode Primer 5	Red	16 µL/tube × 1
	scATAC Barcode Primer 6	Red	16 µL/tube × 1
	scATAC Barcode Primer 7	Red	16 µL/tube × 1
	scATAC Barcode Primer 8	Red	16 µL/tube × 1
	scATAC Barcode Primer 9	Red	16 µL/tube × 1
	scATAC Barcode Primer 10	Red	16 µL/tube × 1
	scATAC Barcode Primer 11	Red	16 µL/tube × 1
	scATAC Barcode Primer 12	Red	16 µL/tube × 1
	scATAC Barcode Primer 13	Red	16 µL/tube × 1
	scATAC Barcode Primer 14	Red	16 µL/tube × 1
	scATAC Barcode Primer 15	Red	16 µL/tube × 1
	scATAC Barcode Primer 16	Red	16 µL/tube × 1

Table 2 DNBelab C Series High-throughput Single-cell ATAC Chip Kit (16 RXN)

Name	Component	Specification	Note
DNBelab C Series High-throughput Single-cell ATAC Chip (Cat. No.: 940-000791-00)	scATAC Chip	16 EA/box × 1	Used together with Single Cell Droplet Generator (Cat. No.: 900-000637-00)
	Sealing Gasket	5 pcs/bag × 1	

1.5 Storage and transportation condition

Table 3 Transportation and storage condition


Name	Storage temperature	Transportation temperature	Validity period
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 1 Droplet Formation Kit) (Cat. No.: 940-000792-00)	2 °C to 8 °C (36 °F to 46 °F)	2 °C to 8 °C (36 °F to 46 °F)	Refer to the label on the box.
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 2 Droplet Formation Kit) (Cat. No.: 940-000794-00)	-25 °C to -15 °C (-13 °F to 5 °F)	-80 °C to -15 °C (-112 °F to 5 °F)	
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 3 Library Preparation Kit) (Cat. No.: 940-000910-00)	-25 °C to -15 °C (-13 °F to 5 °F)	-80 °C to -15 °C (-112 °F to 5 °F)	
DNBelab C Series High-throughput Single-cell ATAC Chip (Cat. No.: 940-000791-00)	10 °C to 30 °C (50 °F to 86 °F)	0 °C to 30 °C (32 °F to 86 °F)	



- Tips**
- The kits are transported by using dry ice. Please check whether dry ice remains upon receipt.
 - When the product is transported, stored, and used appropriately, all of the components retain full activity within the validity period.

1.6 Self-provided materials

Table 4 Self-provided materials

Type	Name	Recommended brand	Cat. No.
Equipment	Single Cell Droplet Generator	MGI	900-000637-00
	Clean bench	/	/
	Fluorescence microscope		
	 Tips A fluorescence microscope is required for counting the number of nuclei.		
	Electronic balance	/	/
	Vortex mixer	/	/
	Mini centrifuge	/	/
	Manual single-channel pipette, with a measurement range as follows: <ul style="list-style-type: none"> • 0.1 µL to 2.5 µL • 0.5 µL to 10 µL • 2 µL to 20 µL • 10 µL to 100 µL • 20 µL to 200 µL • 100 µL to 1000 µL 	/	/
	Manual 8-channel pipette, with a measurement range as follows: <ul style="list-style-type: none"> • 1 µL to 10 µL • 2 µL to 10 µL • 5 µL to 50 µL • 20 µL to 200 µL 	/	/
	Deep-well PCR device (100 µL solution, with heated lid)	/	/
	Centrifuge or equivalents	Eppendorf	5810R
	Magnetic separation rack for 1.5 mL tubes	Thermo Fisher	12321D
	Magnetic separation rack for 0.2 mL tubes	New England Biolabs	S1515S
	Qubit 3.0 fluorometer or equivalents	Thermo Fisher	Q33216

Type	Name	Recommended brand	Cat. No.
Equipment	Agilent 2100 Bioanalyze or equivalents	Agilent Technologies	G2939AA
Reagents	DNA-OFF SOLUTION	TAKARA	9036
	RNase Zap	AMBION	AM9782
	75% ethanol	/	/
	PBS, pH 7.4	Gibco	10010031
	BSA (Bovine Serum Albumin)	Sangon Biotech	A600903-0010
	DAPI or equivalents (Used to dye nuclei)	Sigma-Aldrich	D9542
	Trypan Blue Solution, 0.4% or equivalents	Gibco	15250061
	Nuclease-free water (NF Water)	Ambion	AM9937
	TE buffer, pH 8.0	Ambion	AM9858
	Absolute ethanol (analytical grade)	/	/
	Qubit ssDNA Assay Kit	Invitrogen	Q10212
	Qubit dsDNA HS Assay Kit	Invitrogen	Q32854
	Agilent High Sensitivity DNA Kit or equivalents	Agilent Technologies	5067-4626
Consumables	1 mL syringe	/	/
	0.22 µm filter membrane	PALL	4612
	C-Chip disposable hemocytometer or universal hemocytometer	INCYTO	DHC-N01
	Low-binding sterile filter tips, boxed, with a capacity of 10 µL, 20 µL, 100 µL, 200 µL, or 1000 µL	Axygen	/
	Universal low-binding tips with a capacity of 10 µL, 20 µL, 100 µL, 200 µL, or 1000 µL	Axygen	/
	200 µL wide-bore tips	Axygen	T-205-WB-C
	0.2 mL low-binding PCR tube	Axygen	PCR-02-L-C
	1.5 mL low-binding centrifuge tube	Eppendorf	0030108051
	0.2 mL PCR tube	Axygen	PCR-02-C
	1.5 mL centrifuge tube	Axygen	MCT-150-C

Type	Name	Recommended brand	Cat. No.
Consumables	15 mL centrifuge tube	CORNING	430791
	50 mL centrifuge tube	CORNING	430291
	Qubit Assay Tubes or 0.5mL transparent thin wall PCR tube	Invitrogen/ Axygen	Q32856 or PCR-05-C

1.7 Precautions

- This product is for research use only, and cannot be used for clinical diagnosis. Please read this instructions for use carefully before use.
- Before the experiment, make sure to be familiar with and master the operation methods and precautions of various to-be-used devices.
- Adjust and optimize library preparation procedures according to specific experiment design, sample characteristics, sequencing applications, and devices. The experiment procedures provided in this instructions for use are universal, and the reaction parameters can be adjusted as needed to achieve high performance and efficiency.
- Take out all components of the preparation set in advance, tap to mix the enzymes, briefly centrifuge them, and keep them on ice for further use. Thaw other components on ice, invert them up and down several times after thawing to thoroughly mix them, briefly centrifuge, and place them on ice for further use.
- To avoid an experiment failure caused by cross contamination, it is recommended to perform experiment operations such as sample processing, nuclei transposition, droplet formation, PCR in droplets, demulsification, library preparation, and circularization in a clean laboratory, use low-binding filter tips, and change tips for aspirating different samples.
- It is recommended to proceed with reaction steps in a PCR device with a heated lid mode. The PCR device should be preheated to a required reaction temperature before use.
- Avoid aerosol contamination caused by improper operations on PCR products, which may reduce the accuracy of experimental results. It is recommended to physically separate the PCR reaction solution preparation area from the PCR product purification and detection area. Use special pipettes or other devices, and regularly clean experimental areas (wipe by using 0.5% sodium hypochlorite or 10% bleach) to ensure the cleanliness of the experimental environment.

- Direct contact with skin and eyes should be avoided for all samples and reagents. Do not swallow. Once this happens, immediately rinse with a large amount of water and go to the hospital in time.
- All samples and wastes should be disposed of in accordance with relevant regulations.
- Any other questions, contact the technical support at MGI-service@mgi-tech.com.

Chapter 2 Sample requirements and processing

This chapter describes sample processing and requirements, including precautions and preparation before experiment, sample requirements, experiment requirements, and preparation of reagents and nuclei suspension.

2.1 Precautions before experiment

- It is recommended to operate single-cell experiments in a class 100,000 or 300,000 clean laboratory, or in a clean bench.
- Avoid exogenous nucleic acid contamination when performing single-cell DNA experiment in the clean bench.
- Wear protective equipment such as mask and disposable powder-free latex gloves when performing experiments. During the operation, do not expose the skin of wrist. If the gloves touch the area outside the clean bench, carefully wipe the surface of the gloves with RNase-Zap before continuing the experiment.
- Place all samples on ice in the experiment.
- Consumables such as pipette tips, centrifuge tubes, and sterile water should be sterile, nucleic acid-free, and nuclease-free, and cannot be used for other purposes. Tips should be low-binding nuclease-free filter tips.

2.2 Pre-experiment preparation

2.2.1 Sample requirements

Table 5 Sample requirements

Cell type	Organism/cell samples of human or mammals
-----------	---

Nuclei size	Recommended diameter: < 40 µm
Recommended total nuclei input	<ul style="list-style-type: none"> Recommended total input for transposition : 50,000 to 100,000 nuclei Recommended total input after transposition on chip: 5,000 to 20,000 nuclei
Requirements of nuclei extracted from cell	<ul style="list-style-type: none"> Before nuclei extraction: <ul style="list-style-type: none"> Cytoactivity > 80% Clumping rate < 5% Impurity rate < 5% After nuclei extraction: <ul style="list-style-type: none"> Cytoactivity < 5% Clumping rate < 5% Impurity rate < 5%
Requirements of nuclei extracted from organism	<ul style="list-style-type: none"> Cytoactivity < 5% Clumping rate < 5% Impurity rate < 5%

2.2.2 Experiment requirements

- Before experiment, carefully wipe the gloves, pipettes, bench, and devices with RNase-Zap.
- If a clean bench is used, turn on the light of the clean bench in advance and perform the following steps:
 - Wipe the device and operating deck of the clean bench with DNA-OFF, especially the metal and plastic surfaces.
 - Wait for 10 minutes for degrading the DNA.
 - Turn off the light, and turn on the UV lamp for sterilization for at least 15 minutes.
 - Turn on the light and ventilator after sterilization.

2.2.3 Preparing reagents


Prepare the following reagents:

- PBS (containing 10% BSA)

Table 6 PBS (containing 10% BSA) preparation

Component	Volume
BSA powder	1 g
PBS	Fixed volume to 10 mL

After the BSA powder fully dissolves, use a syringe and 0.22 μ m filter membrane to filter the solution.


 **Tips** The PBS (containing 10% BSA) can be stored at -25 °C to -15 °C (-13 °F to 5 °F) for up to 6 months.

- PBS (containing 1% BSA)

Table 7 PBS (containing 1% BSA) preparation


Component	Volume
PBS	45 mL
PBS (containing 10% BSA)	5 mL

Add the components with the required volumes in proportion and mix them thoroughly.

 **Tips** The PBS (containing 1% BSA) can be stored at 2 °C to 8 °C (36 °F to 46 °F) for up to 1 month.

2.3 Preparing nuclei suspension

For different samples such as a cell-line or live tissue sample, prepare single nuclei suspension in an appropriate way, and wash the nuclei once or twice with the PBS (containing 1% BSA). Resuspend nuclei with an appropriate volume of the PBS (containing 1% BSA) to obtain nuclei suspension. Detect and record the concentration of the nuclei suspension.

-  **Tips**
- It is recommended to mix the nuclei suspension thoroughly by pipetting with a wide-bore tip.
 - When nuclei are counted by using the C-Chip disposable hemocytometer or the universal hemocytometer, a counting result must be correct, or it may affect a final yield. It is recommended to repeat the counting step at least 3 times.
 - After preparation of nuclei suspension is completed, perform the nuclei transposition in 30 minutes.

Chapter 3 Nuclei transposition

This chapter describes how to transpose nuclei by using DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set. The whole procedure takes about 40 minutes.

3.1 Pre-experiment preparation

Table 8 Required reagents

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 2 Droplet Formation Kit) (Cat. No.: 940-000794-00)	Tn5 Transposase	Green	64 μ L/tube \times 1
	5 \times TAG Buffer	Green	80 μ L/tube \times 1


3.2 Nuclei transposition

Perform the following steps:

1. Gently pipette the nuclei suspension prepared in *Preparing nuclei suspension on Page 10* by using a pipette to mix it thoroughly. Take out a 1.5 mL low-binding centrifuge tube, and prepare the transposition reaction solution according to the following table on ice.

Table 9 Transposition reaction solution preparation


Component	Volume (μ L) required for each tube
5 \times TAG Buffer	5
Tn5 Transposase	4
PBS (containing 1% BSA)	16-X
Nuclei suspension	X
Total	25

-  **Tips**
- The recommended total number of input nuclei ranges from 50,000 to 100,000, and the recommended concentration of the nuclei is 6,250 nuclei/ μL $< N < 10,000$ nuclei/ μL . The biggest volume of the nuclei should be 16 μL and may vary according to the nuclei concentration. If the volume is less than 16 μL , the PBS (containing 1% BSA) is added as a complement.
 - Tn5 Transposase should be tapped gently and centrifuged to the bottom of the tube briefly before use. If not mixed thoroughly, it may cause in failure of the experiment.
 - Ensure that all reagents are added to the bottom of the centrifuge tube. Failure to do so may cause the reagents to remain on the tube wall and cause in loss of reaction solution.

2. After the reaction solution is prepared, adjust the pipette range to 20 μL . Gently pipette to mix the reaction solution thoroughly.

-  **Tips** Gently pipette and tap the tube, ensure that no bubbles exist in the tube.

3. Put the centrifuge tube in the mixer at 37 $^{\circ}\text{C}$ (99 $^{\circ}\text{F}$), and transpose at 500 rpm for 30 minutes. After transposing for 15 minutes, take out the tube from the mixer, pipette to mix it thoroughly, and then place the tube in the mixer to continue transposing.
4. After the transposition is completed, place the 1.5 mL centrifuge tube on ice immediately and gently pipette to mix it thoroughly. Count and records the number of the nuclei.

-  **Tips** After the transposition is completed, perform the droplet formation in 30 minutes.

Chapter 4 Droplet formation

This chapter describes how to form droplets from the prepared nuclei suspension after transposition by using DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set. The whole procedure takes about 20 minutes.

4.1 Pre-experiment preparation

4.1.1 Preparing reagents and consumables

Table 10 Required reagents and consumables

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 1 Droplet Formation Kit) (Cat. No.: 940-000792-00)	scATAC Cell Beads	Black	1.6 mL/tube × 1
	scATAC Wash Buffer	Brown	1.6 mL/tube × 2
	scATAC Lysis Buffer	Black	768 µL/tube × 1
	P100 Oil	Natural	11.2 mL/bottle × 1
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 2 Droplet Formation Kit) (Cat. No.: 940-000794-00)	scATAC Bead Buffer	Black	1.114 mL/tube × 1
	PCR Primer	Black	39 µL/tube × 1
	PCR Enzyme I	White	128 µL/tube × 1
	scATAC Cell Buffer	White	748 µL/tube × 1
	scATAC DIR Reagent	White	32 µL/tube × 1
DNBelab C Series High-throughput Single-cell ATAC Chip (Cat. No.: 940-000791-00)	scATAC Chip	/	16 EA/box* 1
	Sealing Gasket	/	5 pcs/bag* 1



- Tips**
- Take out the P100 Oil for at least 30 minutes in advance to equilibrate to room temperature.
 - Take out the scATAC Lysis Buffer in advance to equilibrate to room temperature until all the crystals in the solution are dissolved.
 - DNBELAB C Series High-throughput Single-cell ATAC Chip is used together with Single Cell Droplet Generator (Cat. No.: 900-000637-00).

4.1.2 Preparing sample phase suspension

Perform the following steps:

1. Gently pipette the nuclei suspension prepared in *Nuclei transposition on Page 11* by using a pipette to mix it thoroughly, and prepare the sample phase suspension according to the following table.

Table 11 Sample phase suspension preparation

Component	Volume (μL) required for each tube
scATAC Cell Buffer	46.7
PCR Enzyme I	8
scATAC DIR Reagent	2
NF Water	43.3-X
Nuclei suspension	X
Total	100



Tips • 'X' indicates the volume of the nuclei suspension.


- The recommended total number of input nuclei suspension ranges from 5,000 to 20,000, and the biggest volume of the nuclei suspension should be 43.3 μL and may vary according to the nuclei concentration. If the volume is less than 43.3 μL, the NF Water is added as a complement.
 - If there are identical samples that need to use multiple chips, these samples can be combined to make the sample phase suspension.
 - You can add the lost reagents in the process when making the sample phase suspension.
 - For the resuscitated or fragile nuclei, it is recommended to use a wide-bore tip to pipette and mix the suspension thoroughly.
2. After the suspension is prepared, adjust the pipette range to 80 μL, gently pipette the mixture to mix it thoroughly, and place it on ice for future use.

4.1.3 Preparing beads phase suspension

Perform the following steps:

1. Take out the scATAC Cell Beads, and gently invert or pipette it to mix it thoroughly.
2. Aspirate 100 μL of the scATAC Cell Beads (for each sample), and transfer it to a 0.2 mL low-binding centrifuge tube (the volume varies according to the number of samples).
3. Place the PCR tube on the magnetic separation rack for 3 to 5 minutes. Gently remove and dispose of the supernatant to avoid loss of the beads.
4. Add 200 μL of the scATAC Wash Buffer (for each sample) and clean the beads by rotating the centrifuge tube 180 degrees on the magnetic separation rack. After the supernatant becomes clear, use a pipette to remove and dispose of the supernatant.

5. Remove the PCR tube from the magnetic separation rack, and add 69.6 μL of the scATAC Bead Buffer, 48 μL of the scATAC Lysis buffer, 2.4 μL of the PCR Primer accordingly to prepare the beads phase suspension.
6. Adjust the pipette range to 100 μL , gently pipette the mixture to mix it thoroughly, and place it on ice for future use.

 **Tips** If there are multiple samples that need to use multiple chips, it is recommended to prepare the beads phase suspension required for all the samples in a 1.5 mL low-binding centrifuge tube.

4.2 Performing droplet formation

Perform the following steps:

1. Take out the scATAC Chip and the single cell droplet generator.




Figure 1 scATAC chip




Figure 2 Single cell droplet generator (left: front view, right: back view)

2. Switch the power switch of the single cell droplet generator to **Open** position. Tap **Open**, the compartment door opens automatically. Take out the chip stage.

3. Unlock the buckle of the chip stage, install the chip by pushing it into the chip stage from the side (the number of the loaded chip depends on the number of the sample), and lock the buckle at the oil well position.




Figure 3 Loading the chip

4. Snap the sealing gasket into the collection well (tailor the gasket according to the number of the chip).




Figure 4 Loading the sealing gasket

5. Place the loaded chip stage into the compartment.




Figure 5 Loading the loaded chip stage

6. Add 100 μL of the sample phase suspension, 700 μL of the P100 oil, 120 μL of the beads phase suspension in order accordingly to the sample well, oil well and beads well.



- Tips**
- Before adding the sample phase suspension and beads phase suspension, separately pipettes them to mix thoroughly. Ensure that no bubbles exist during pipetting. When adding the solutions, do not suspend the tip. Gently add the sample near the edge of the well to avoid forming bubbles at the bottom.
 - The total time for adding the three solutions should be within 1 minute.
 - Ensure that you add the sample phase suspension first, then add the P100 Oil, and add the beads phase suspension at last. Strictly follow this order. Otherwise, it may cause in failure of droplet formation.




Figure 6 Adding solutions to the chip

7. Tap **Close** to close the compartment door.
8. Select **ATAC** as the reaction type, and select the corresponding channel.



- Tips** The channel **A/B/C/D** on the screen corresponds to the channel "A"/"B"/"C"/"D" on the stage.

9. Tap **Start** to start the reaction.




Figure 7 Droplet formation starts

10. After the reaction is completed, tap **Confirm**.





Figure 8 Droplet formation completed

11. Tap **Open** to open the compartment door and take out the chip stage. Remove the sealing gasket, unlock the buckle, take out the chip, and recover the droplets.

 **Tips** After the droplet formation is completed, recover the droplets immediately to avoid droplets evaporation due to exposure to the air or loss of droplets due to attachment on the well.

Chapter 5 Performing PCR in droplets

This chapter describes how to recover droplets after the droplet formation is completed and how to perform PCR in droplets. The whole procedure takes about 1 hour and 20 minutes.

5.1 Pre-experiment preparation

Table 12 Required reagents

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 1 Droplet Formation Kit) (Cat. No.: 940-000792-00)	Mineral Oil	Brown	6.4 mL/bottle × 1



Tips Take out the Mineral Oil at least 30 minutes in advance and equilibrate to room temperature.

5.2 Recovering droplets

Perform the following steps:

1. Take out a clean PCR 8-strip tube and mark it with a marker.



Tips Use 1 PCR 8-strip tube for 1 sample.

2. Gently aspirate all the droplets in the collection well of the chip by using a 200 μ L low-binding tip (you can slightly lean the chip to collect the droplets). Suspend the pipette vertically for a few seconds until the droplets float to the upper level of the oil. Gently transfer the oil at the bottom to the last 4 tubes of the PCR 8-strip tube and the droplets to the first 4 tubes of the PCR 8-strip tube. Repeat this step until all the droplets have been transferred.



- Tips**
- The first 4 tubes of the PCR 8-strip tube are loaded with all the droplets of 1 sample, and the last 4 tubes are loaded with P100 Oil. Transfer the droplets gently. Pipetting or tapping the droplets intensively will break the droplets.
 - Divide the droplets of 1 sample into the first 4 tubes of the PCR 8-strip tube on average. The volume of the droplets in each tube should be about 50 μ L to 100 μ L.

3. Use the oil collected in the last 4 tubes of the PCR 8-strip tube to rinse the remained droplets in the collection well (do not pipette). Transfer the droplets to the first 4 tubes of the PCR 8-strip tube. Try to recover all the droplets.



Tips Do not place the droplets after formation for over 40 minutes. Otherwise, the data quality may be influenced.

4. Add 100 μ L of the Mineral Oil on the top level of the first 4 tubes of the PCR 8-strip PCR tube.




Figure 9 Liquid in the 8-strip tube

Tips The liquid in the PCR 8-strip tube is divided into 3 levels: the upper level is the mineral oil, the middle level is the droplets, and the bottom level is P100 Oil (little or none).

5. Perform PCR in droplets in the PCR device according to the following table.

Table 13 PCR in droplets condition (100 μ L reaction solution)

Temperature	Time	Cycles
105 °C (221 °F) (heated lid)	On	/
72 °C (162 °F)	5 min	1
98 °C (208 °F)	30 s	1
98 °C (208 °F)	10 s	10
63 °C (145 °F)	30 s	
72 °C (162 °F)	1 min	
72 °C (162 °F)	5 min	1
4 °C (39 °F)	Hold	/

Stop point After the PCR reaction is completed, the solution can store at 2 °C to 8 °C (36 °F to 46 °F) for 72 hours.


Chapter 6 Demulsification and beads selection

This chapter describes the procedure of recovering DNA by demulsification. The whole procedure takes about 1 hour.

6.1 Pre-experiment preparation

Table 14 Required reagents


Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 1 Droplet Formation Kit) (Cat. No.: 940-000792-00)	scATAC Breakage Reagent	Brown	1.6 mL/tube × 1
	DNA Clean Beads	Natural	5.933 mL/bottle × 1

-  **Tips**
- Take out the DNA Clean Beads at least 30 minutes in advance to equilibrate to room temperature. Before use, mix it thoroughly.
 - Use the scATAC Breakage Reagent in the ventilation cabinet.


6.2 Performing demulsification

Perform the following steps:


1. Transfer the droplets in the middle level of the first 4 tubes of the PCR 8-strip tube prepared in *Recovering droplets on Page 19* to a clean 1.5 mL low-binding centrifuge tube.

 **Tips** Avoid pipetting the upper level of the Mineral Oil. Pipetting the P100 Oil in the bottom level will not make an affect.

2. Add 100 μ L of the scATAC Breakage Reagent in the centrifuge tube, invert the centrifuge tube 15 to 20 times, and place it at room temperature for 3 minutes.

 **Tips** Do not shake or vortex the tube strongly.

3. Centrifuge the centrifuge tube for 1 minute at room temperature with a centrifugal force of 1000 \times g.
4. Place the centrifuge tube on the magnetic separation rack for 3 to 5 minutes.
5. Slowly pipette 100 μ L of the aqueous phase to a new 1.5 mL low-binding centrifuge tube.

 **Tips** Slowly pipette the aqueous phase. Avoid pipetting the oil phase in the bottom level, the interface layer between the aqueous phase and the oil phase in the middle and the remained mineral oil in the upper level.

6.3 Beads selection



- Tips**
- Take out the DNA Clean Beads at least 30 minutes in advance to equilibrate to room temperature. Before use, mix it thoroughly.
 - Before operation, carefully read through *Appendix 1 About the DNA Clean Beads and purification on Page 36*.

Perform the following steps:

1. Aspirate 50 μ L of the DNA Clean Beads into the low-binding centrifuge tube with the aqueous phase prepared in *Performing demulsification on Page 21* and pipette to mix it thoroughly. In the last pipetting, ensure that all the liquids and beads in the pipette are transferred into the centrifuge tube.
2. Incubate the centrifuge tube at room temperature for 5 minutes.
3. After centrifuge briefly, place and keep the centrifuge tube on the magnetic separation rack for 2 to 5 minutes until the liquid becomes clear. Gently transfer the supernatant to a new 1.5 mL low-binding centrifuge tube by using a pipette.



- Tips**
- Do not dispose of but reserve the supernatant in this step.
 - Do not pipette the beads.
4. Aspirate 50 μ L of the DNA Clean Beads to the centrifuge tube containing the supernatant and gently pipette and mix it thoroughly. In the last pipetting, ensure that all the liquid and beads in the pipette are transferred into the PCR tube.
 5. Incubate the PCR tube at room temperature for 5 minutes.
 6. After centrifuge briefly, place and keep the centrifuge tube on the magnetic separation rack for 2 to 5 minutes until the liquid becomes clear. Gently remove and dispose of the supernatant.
 7. Keep the centrifuge tube on the magnetic separation rack, add 500 μ L of the 80% ethanol freshly prepared, rinse the beads and the tube wall, and keep the tube for 30 seconds. Gently remove and dispose of the supernatant.
 8. Repeat the last step, try to remove all the liquid from the tube. If some residuals remain on the tube wall, briefly centrifuge the tube, place it on the magnetic separation rack for separation, and remove all liquid at the bottom of the tube by using a pipette with a small measurement range.
 9. Keep the tube on the magnetic rack, and open the cap of the tube to dry the beads at room temperature until no wetness or glossiness is visible on the beads surface. Do not let the beads overdried to avoid beads from cracking.
 10. Remove the centrifuge tube from the magnetic separation rack, add 22 μ L of the NF water to elute DNA, and pipette to mix it thoroughly.
 11. Incubate the centrifuge tube at room temperature for 5 minutes.

12. Centrifuge the tube briefly, place the tube on the magnetic separation rack, and wait for 2 to 5 minutes until the liquid becomes clear. Transfer 20 μ L of the supernatant to a new 1.5 mL centrifuge tube and mark it as 'Product'.
13. Take out 1 μ L of the 'Product' and detect its concentration of product by using the Qubit dsDNA HS Assay Kit.

II Stop point The 'Product' can be stored at -25 °C to -15 °C (-13 °F to 5 °F) for 6 months.

Chapter 7 Library preparation


This chapter describes the library preparation. The whole procedure takes about 1 hour and 20 minutes.

7.1 Pre-experiment preparation

Table 15 Required reagents

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 1 Droplet Formation Kit) (Cat. No.: 940-000792-00)	DNA Clean Beads	Natural	5.933 mL/bottle × 1
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 3 Library Preparation Kit)) (Cat. No.: 940-000910-00)	PCR Enzyme II	Blue	800 μ L/tube × 1
	scATAC Barcode Primer 1	Red	16 μ L/tube × 1
	scATAC Barcode Primer 2	Red	16 μ L/tube × 1
	scATAC Barcode Primer 3	Red	16 μ L/tube × 1
	scATAC Barcode Primer 4	Red	16 μ L/tube × 1
	scATAC Barcode Primer 5	Red	16 μ L/tube × 1
	scATAC Barcode Primer 6	Red	16 μ L/tube × 1
	scATAC Barcode Primer 7	Red	16 μ L/tube × 1
	scATAC Barcode Primer 8	Red	16 μ L/tube × 1
	scATAC Barcode Primer 9	Red	16 μ L/tube × 1

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 3 Library Preparation Kit)) (Cat. No.: 940-000910-00)	scATAC Barcode Primer 10	Red	16 µL/tube × 1
	scATAC Barcode Primer 11	Red	16 µL/tube × 1
	scATAC Barcode Primer 12	Red	16 µL/tube × 1
	scATAC Barcode Primer 13	Red	16 µL/tube × 1
	scATAC Barcode Primer 14	Red	16 µL/tube × 1
	scATAC Barcode Primer 15	Red	16 µL/tube × 1
	scATAC Barcode Primer 16	Red	16 µL/tube × 1

 **Tips** Take out the DNA Clean Beads at least 30 minutes in advance to equilibrate to room temperature.

7.2 Performing library preparation

 **Tips** Please read through *Appendix 2 Using scATAC Barcode Primer (1 to 16) on Page 37*.

Perform the following steps:

1. According to the 'Product' concentration detected in *Demulsification and beads selection on Page 20*, add 3 ng (if less than 3 ng, put them all in) of the 'Product' in a new 0.2 mL PCR tube, add NF Water as a complement to a total volume of 46 µL, and place the PCR tube on ice.
2. Add 4 µL of the scATAC Barcode Primer in the PCR tube.

 **Tips** Record the number of scATAC Barcode Primer added in each sample.

3. Add 50 µL of the PCR Enzyme II, vortex to mix it thoroughly, and centrifuge the tube briefly to collect the reaction solution to the bottom of the tube.
4. Place the PCR tube on the PCR device and perform the reaction according to the following table.

Table 16 PCR reaction condition (100 μ L reaction solution)

Temperature	Time	Cycles
105 °C (221 °F) (heated lid)	On	/
98 °C (208 °F)	30 s	1
98 °C (208 °F)	10 s	12
63 °C (145 °F)	30 s	
72 °C (162 °F)	1 min	
72 °C (162 °F)	5 min	1
4 °C (39 °F)	Hold	/

7.3 Beads selection



Tips Before operation, carefully read through *Appendix 1 About the DNA Clean Beads and purification on Page 36*.

Perform the following steps:

1. Aspirate 40 μ L of the DNA Clean Beads into the centrifuge tube with the PCR products prepared in *Performing library preparation on Page 24* and gently pipette and mix it thoroughly. In the last pipetting, ensure that all the liquids and beads in the pipette are transferred into the PCR tube.
2. Incubate the PCR tube at room temperature for 5 minutes.
3. After centrifuge briefly, place and keep the PCR tube on the magnetic separation rack for 2 to 5 minutes until the liquid becomes clear. Gently transfer the supernatant to a new 1.5 mL low-binding centrifuge tube by using a pipette.



Tips

- Do not dispose of but reserve the supernatant in this step.
- Do not pipette the beads.

4. Aspirate 70 μ L of the DNA Clean Beads to the centrifuge tube containing the supernatant prepared in the last step and gently pipette and mix it thoroughly. In the last pipetting, ensure that all the liquid and beads in the pipette are transferred into the centrifuge tube.
5. Incubate the centrifuge tube at room temperature for 5 minutes.
6. After centrifuge briefly, place and keep the centrifuge tube on the magnetic separation rack for 2 to 5 minutes until the liquid becomes clear. Gently remove and dispose of the supernatant.

7. Place and keep the centrifuge tube on the magnetic separation rack, add 500 μL of the 80 % ethanol freshly prepared to rinse the beads and the tube wall, keep the tube for 30 seconds, and dispose of the supernatant.
8. Repeat the last step, try to remove all liquid from the tube. If some residuals remain on the tube wall, briefly centrifuge the tube, place it on the magnetic separation rack for separation, and remove all liquid at the bottom of the tube by using a pipette with a small measurement range.
9. Keep the tube on the magnetic rack, and uncap the tube to dry the beads at room temperature until no wetness or glossiness is visible on the beads surface. Do not let the beads overdried to avoid the beads from cracking.
10. Remove the centrifuge tube from the magnetic separation rack, add into 42 μL of the TE Buffer, and pipette to mix it thoroughly.
11. Incubate the centrifuge tube at room temperature for 5 minutes.
12. Centrifuge the tube briefly, place the tube on the magnetic separation rack, and wait for 2 to 5 minutes until the liquid becomes clear. Transfer 40 μL of the supernatant to a new 1.5 mL centrifuge tube and mark it as 'Library'.
13. Take out 1 μL of the 'Library' and detect the concentration by using the Qubit dsDNA HS Assay Kit. Take out 1 μL of the 'Library' and detect the segment distribution by using the Agilent High Sensitivity DNA Kit. The following figure is an example of the library segment distribution:






Figure 10 Library segment distribution (Agilent 2100)

-  **Tips** Reference value: When the 'Product' input is 3 ng, library concentration should be greater than 50 ng/ μL .
-  **Stop point** The 'Library' can be stored at -25 $^{\circ}\text{C}$ to -15 $^{\circ}\text{C}$ (-13 $^{\circ}\text{F}$ to 5 $^{\circ}\text{F}$) for 6 months.

Chapter 8 Circularization

This chapter describes the whole procedure of single-stranded circularization and products purification, which takes about 1 hour and 30 minutes.

8.1 Pre-experiment preparation

Table 17 Required reagents

Name	Component	Cap color	Specification
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 1 Droplet Formation Kit) (Cat. No.: 940-000792-00)	DNA Clean Beads	Natural	5.933 mL/bottle × 1
DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 3 Library Preparation Kit) (Cat. No.: 940-000910-00)	scATAC Splint Buffer	Yellow	80 µL/tube × 1
	scATAC Cyclization Buffer	Yellow	157 µL/tube × 1
	scATAC DNA Ligase	Yellow	8 µL/tube × 1
	scATAC Exo Enzyme	Orange	42 µL/tube × 1
	scATAC Exo Buffer	Orange	23 µL/tube × 1
	scATAC Stop Buffer	Orange	48 µL/tube × 1



Tips Take out the DNA Clean Beads at least 30 minutes in advance to equilibrate to room temperature. Vortex to mix it thoroughly before use.

8.2 Denaturation

Perform the following steps:

1. Add 400 ng of the 'Library' into a new 0.2 mL PCR tube and add TE Buffer as a complement to a total volume of 45 µL.
2. Add 5 µL of the scATAC Splint Buffer into the PCR tube, vortex to mix it thoroughly, and centrifuge briefly to collect the reaction solution to the bottom of the tube.
3. Put the PCR tube on the PCR device and perform the reaction according to the following table.

Table 18 Denaturation condition (50 µL reaction solution)

Temperature	Time
105 °C (221 °F) (heated lid)	On
95 °C (203 °F)	3 min

- After the reaction is completed, place and keep the PCR tube on ice for 5 minutes, and then add single-stranded circularization reaction solution into the tube.



- Tips**
- After the reaction at 95 °C (203 °F) is completed, place the PCR tube on ice immediately. Otherwise, circularization efficiency may be affected by the DNA renaturation.
 - You can prepare the single-stranded circularization reaction solution during the waiting period of 5 minutes.

8.3 Single-stranded circularization

Perform the following steps:

- Prepare the reaction solution for single-stranded circularization on ice according to the following table.

Table 19 Single-stranded circularization reaction solution

Component	Volume (µL) required for each tube
scATAC Cyclization Buffer	9.5
scATAC DNA Ligase	0.5
Total	10

- Aspirate 10 µL of the single-stranded circularization reaction solution into the PCR tube, vortex to mix it thoroughly, and centrifuge briefly to collect the reaction solution to the bottom of the tube.
- Put the PCR tube on the PCR device and perform the reaction according to following table.

Table 20 Single-stranded circularization condition (60 µL reaction solution)

Temperature	Time
50 °C (122 °F) (heated lid)	On
37 °C (99 °F)	30 min
4 °C (39 °F)	Hold

4. After the reaction is completed, centrifuge the PCR tube briefly and place the tube on ice. Perform the enzyme digestion immediately.

8.4 Enzyme digestion

Perform the following steps:

1. Prepare the enzyme digestion reaction solution on ice according to the following table.

Table 21 Enzyme digestion reaction solution

Component	Volume (μL) required for each tube
scATAC Exo Enzyme	2.6
scATAC Exo Buffer	1.4
Total	4.0

2. Aspirate the 4 μL enzyme digestion reaction solution into the PCR tube prepared in *Single-stranded circularization on Page 28* by using the pipette, vortex to mix it thoroughly, and centrifuge briefly to collect the reaction solution to the bottom of the tube.
3. Place the PCR tube on the PCR device and perform the reaction according to the following table.

Table 22 Digestion condition (64 μL reaction solution)

Temperature	Time
50 °C (122 °F) (heated lid)	On
37 °C (99 °F)	30 min
4 °C (39 °F)	Hold

4. After the reaction is completed, add 3 μL of the scATAC Stop Buffer into the PCR tube, vortex to mix it thoroughly, and centrifuge it briefly to collect the reaction solution to the bottom of the tube.

8.5 Beads purification



- Tips**
- Take out the DNA Clean Beads at least 30 minutes in advance to equilibrate to room temperature. Vortex to mix it thoroughly before use.
 - Before operation, carefully read through *Appendix 1 About the DNA Clean Beads and purification on Page 36*.

Perform the following steps:

1. Transfer the digestion products prepared in *Enzyme digestion on Page 29* into a new 1.5 mL low-binding centrifuge tube. Add 160.8 μL of the DNA Clean Beads into the centrifuge tube and gently pipette and mix it thoroughly. In the last pipetting, ensure that all the liquids and beads in the pipette are transferred into the centrifuge tube.
2. Incubate the centrifuge tube at room temperature for 10 minutes.
3. Centrifuge briefly, place and keep the centrifuge tube on the magnetic separation rack for 2 to 5 minutes until the liquid becomes clear. Gently remove and dispose of the supernatant.
4. Place and keep the centrifuge tube on the magnetic separation rack. Add 500 μL of the 80 % ethanol freshly prepared to rinse the beads and the tube wall, wait for 30 seconds, and dispose of the supernatant.
5. Repeat the last step, try to remove all liquids from the tube. If some residuals remain on the tube wall, briefly centrifuge the tube, place it on the magnetic separation rack for separation, and remove all liquid at the bottom of the tube by using a pipette with a small measurement range.
6. Keep the tube on the magnetic rack, and uncap the tube to dry the beads at room temperature until no wetness or glossiness is visible on the beads surface. Do not let the beads overdried to avoid the beads from cracking.
7. Remove the centrifuge tube from the magnetic separation rack, add into 32 μL of the TE Buffer, and pipette to mix it thoroughly.
8. Incubate the centrifuge tube at room temperature for 10 minutes.
9. After centrifuge briefly, place and keep the centrifuge tube on the magnetic separation rack for 2 to 5 minutes until the liquid becomes clear. Transfer 30 μL of the supernatant to a new 1.5 mL centrifuge tube gently and mark the tube as 'ssDNA'.



Tips Do not pipette the beads in this step. Failure to do so will take a great effect on the quality of the library sequencing.

10. Take out 1 μL of the 'ssDNA' and detect the concentration by using the Qubit ssDNA Assay Kit.



Tips Reference value: The 'ssDNA' concentration should be greater than 0.5 ng/ μL .



Stop point The 'ssDNA' can be stored at -25 $^{\circ}\text{C}$ to -15 $^{\circ}\text{C}$ (-13 $^{\circ}\text{F}$ to 5 $^{\circ}\text{F}$) for 1 month.

Chapter 9 Sequencing

This chapter describes the sequencer, sequencing reagents and read length that are compatible to the library. The library structure is as shown below:




Figure 11 Library structure

- Tips**
- Read 1 = 115 bp (Read 1 fixed sequence 6 + 6 + 33 = 45 bp dark reaction)
 - Read 2 = 69 bp (Read 2 fixed sequence 19 bp dark reaction)
 - Sample Barcode = 10 bp

9.1 Sequencing requirements of DNBSEQ-G400RS

9.1.1 Pre-experiment preparation

- Tips** According to the sequencing set, carefully read *DNBSEQ-G400RS HotMPS High-throughput Sequencing Set User Manual* or *DNBSEQ-G400RS High-throughput (Rapid) Sequencing Set User Manual* before starting sequencing, and strictly perform all operations according to the operations in the user manual.

Table 23 Preparation list

Type	Model No.	Cat. No.
Genetic sequencer	DNBSEQ-G400RS	/
Sequencing set	HotMPS High-throughput Sequencing Set (G400 HM FCL PE100)	940-000489-00
	HotMPS High-throughput Sequencing Set (G400 HM FCL PE150)	940-000244-00
	DNBSEQ-G400RS High-throughput Sequencing Set (FCL PE100)	1000016950
	DNBSEQ-G400RS High-throughput Sequencing Set (FCL PE150)	1000016952

9.1.2 Making DNB

Make DNB according to the following table.

Table 24 DNBSEQ-G400RS making DNB requirements

Model No.	DNBSEQ-G400RS (HotMPS)		DNBSEQ-G400RS (StandardMPS)	
Sequencing read	PE100	PE150	PE100	PE150
Make DNB input	10 ng	10 ng	10 ng	10 ng
RCA time	20 min	20 min	20 min	20 min

9.1.3 Pooling libraries

For details about protocols for pooling different libraries based on the scATAC Barcode Primer (1 to 16), refer to *Appendix 2 Using scATAC Barcode Primer (1 to 16) on Page 37*.

9.1.4 Sequencing parameters

Table 25 DNBSEQ-G400RS (HotMPS) sequencing software version and read length (mixed sample, sequencing barcode)

Software version	ECR 1.0
Control software version	1.7.1.1799
Basecall version	1.5.0.323
Sequencing script	Z_scATAC_BC_HM
Read 1	115 cycles (1 to 6 bp, 17 to 22 bp, 33 to 65 bp are set for dark reaction)
Read 2	69 cycles (1 to 19 bp are set for dark reaction)
Sample barcode	10 cycles

Table 26 DNBSEQ-G400RS (HotMPS) sequencing software version and read length (unmixed sample, not sequencing barcode)

Software version	ECR 1.0
Control software version	1.7.1.1799
Basecall version	1.5.0.323
Sequencing script	Z_scATAC_noBC_HM


Read 1	115 cycles (1 to 6 bp, 17 to 22 bp, 33 to 65 bp are set for dark reaction)
Read 2	69 cycles (1 to 19 bp are set for dark reaction)
Sample barcode	/

Table 27 DNBSEQ-G400RS (StandardMPS) sequencing software version and read length (mixed sample, sequencing barcode)

Software version	ECR 3.0 or later version
Control software version	DNBSEQ-G400_1.0.0.34 or later version
Basecall version	Basecall_1.0.8.208 or later version
Sequencing script	scATAC_BC
Read 1	115 cycles (1 to 6 bp, 17 to 22 bp, 33 to 65 bp are set for dark reaction)
Read 2	69 cycles (1 to 19 bp are set for dark reaction)
Sample barcode	10 cycles

Table 28 DNBSEQ-G400RS (StandardMPS) sequencing software version and read length (unmixed sample, not sequencing barcode)

Software version	ECR 3.0 or later version
Control software version	DNBSEQ-G400_1.0.0.34 or later version
Basecall version	Basecall_1.0.8.208 or later version
Sequencing script	scATAC_noBC
Read 1	115 cycles (1 to 6 bp, 17 to 22 bp, 33 to 65 bp are set for dark reaction)
Read 2	69 cycles (1 to 19 bp are set for dark reaction)
Sample barcode	/

 **Tips** If the number of nuclei input is 10,000, the recommended sequencing data is 400 M Reads. The average sequencing depth for each nucleus is 50 K and can be modified according to requirements.

9.2 Sequencing requirements of DNBSEQ-T7RS

9.2.1 Pre-experiment preparation


 **Tips** According to the sequencing set, carefully read *DNBSEQ-T7RS HotMPS High-throughput Sequencing Set User Manual* or *DNBSEQ-T7RS High-throughput Sequencing Set User Manual* before starting sequencing, and strictly perform all operations according to the operations in the user manual.

Table 29 Preparation list

Type	Model No.	Cat. No.
Genetic sequencer	DNBSEQ-T7RS	/
Sequencing set	DNBSEQ-T7RS HotMPS High-throughput Sequencing Set (FCL PE100)	940-000247-00
	DNBSEQ-T7RS HotMPS High-throughput Sequencing Set (FCL PE150)	940-000248-00
	DNBSEQ-T7RS High-throughput Sequencing Set (FCL PE100)	1000028455
	DNBSEQ-T7RS High-throughput Sequencing Set (FCL PE150)	1000028454

9.2.2 Making DNB

Make DNB according to the following table.

Table 30 DNBSEQ-T7RS making DNB requirements

Model No.	DNBSEQ-T7RS (HotMPS)		DNBSEQ-T7RS (StandardMPS)	
Sequencing read	PE100	PE150	PE100	PE150
Make DNB input	15 ng	15 ng	10 ng	10 ng
RCA time	10 min	10 min	30 min	15 min

9.2.3 Pooling libraries

For details about protocols for pooling different libraries based on the scATAC Barcode Primer (1 to 16), refer to *Appendix 2 Using scATAC Barcode Primer (1 to 16) on Page 37*.

9.2.4 Sequencing parameters


Table 31 DNBSEQ-T7RS (HotMPS) sequencing software version and read length

Software version	ECR 4.0 or later version
Control software version	1.4.1.812 or later version
Basecall version	BCS_1.0.10.142 or later version
Sequencing script	Customize
Custom primers	No
Read 1	115 cycles (1 to 6 bp, 17 to 22 bp, 33 to 65 bp are set for dark reaction)
Read 2	69 cycles (1 to 19 bp are set for dark reaction)
Sample barcode	10 cycles

Table 32 DNBSEQ-T7RS (StandardMPS) sequencing software version and read length

Software version	ECR 3.0 or later version
Control software version	1.3.3.553 or later version
Basecall version	1.4.2.47_Ubuntu or later version

Sequencing script	Customize
Custom primers	No
Read 1	115 cycles (1 to 6 bp, 17 to 22 bp, 33 to 65 bp are set for dark reaction)
Read 2	69 cycles (1 to 19 bp are set for dark reaction)
Sample barcode	10 cycles

 **Tips** If the number of nuclei input is 10,000, the recommended sequencing data is 400 M Reads. The average sequencing depth for each nucleus is 50 K and can be modified according to requirements.

Appendix 1 About the DNA Clean Beads and purification

Precautions before use of the DNA Clean Beads

- To ensure recovery efficiency of the DNA Clean Beads (hereinafter called the beads), take it out from the 4 °C (39 °F) refrigerator 30 minutes in advance, vortex it to mix thoroughly and equilibrate to room temperature.
- Vortex or invert the beads before use each time, to ensure that the beads are mixed thoroughly.
- The volume of the beads directly affects the lower limit length of the purified DNA fragments. The higher volume of the beads leads to the smaller lower limit length of the purified DNA fragments.

Precautions during use of the DNA Clean Beads

- If the volume of to be purified DNA segments decreases due to evaporation caused by incubation, the TE Buffer should be added to complement the liquid to the required volume, and then the recommended volume of the beads is added to perform purification.
- After mixing the DNA segments and the beads thoroughly and placing the centrifuge tube on the magnetic separation rack for separation, remove the supernatant after the solution becomes completely clear. This process usually

takes 2 to 5 minutes. Because magnetism of magnetic separation racks might be different, the separation time may be longer, depending on the time in which the solution becomes completely clear.

- When separating the beads from the solution, avoid contact between the pipette tip and the beads. Reserve 2 or 3 μL of the solution, to avoid aspirating the beads. If the beads are aspirated accidentally, pipette down all the beads and the liquid into the centrifuge tube, and re-aspirate the supernatant after separation.
- Rinse the beads by using 80% ethanol that is freshly prepared and equilibrated to room temperature. During rinsing, keep the centrifuge tube on the magnetic separation rack, and operate the pipette tip on the side that is away from the magnetic separation rack. Do not pipette and stir the beads.
- Try to remove all liquid from the tube in second rinsing with the ethanol. If some residuals remain on the tube wall, briefly centrifuge the tube, place it on the magnetic separation rack for separation, and remove all liquid at the bottom of the tube by using a pipette with a small measurement range.
- After rinsing the beads with the ethanol twice, completely dry the beads at room temperature. Incomplete drying of the beads (with a reflective surface) easily causes absolute ethanol residuals and affects subsequent reactions. Over-drying of the beads (beads with cracks) reduces the purification yield. Drying at room temperature usually takes 5 to 10 minutes. The drying time varies with the room temperature and humidity. Start elution after observing that the beads are not reflective.
- Avoid contact between the pipette tip and the beads when removing the supernatant after elution, as removal of the beads might affect subsequent purification reactions. Therefore, the elution volume should be 2 μL greater than the volume of the aspirated supernatant.
- Carefully open or close the lid of the 1.5 mL centrifuge tube on the magnetic separation rack, to avoid spill of the beads or liquid caused by strong shaking. It is recommended to hold the tube at the middle and lower part with your fingers when opening the lid.

Appendix 2 Using scATAC Barcode Primer (1 to 16)

DNBelab C Series High-throughput Single-cell ATAC Library Preparation Set (Box 3 Library Preparation Kit) includes 16 tubes of scATAC Barcode Primer. The scATAC Barcode Primers in the preparation set are designed based on the base balancing principle. To ensure high performance, carefully read through this appendix.



- Tips**
- Avoid placing the tubes at a temperature higher than the room temperature. Otherwise, melting occurs, which affects the performance.
 - Mix each tube of scATAC Barcode Primer thoroughly and centrifuge it before use, and wipe the lid of the tube with the lint-free paper. Gently open the lid during use, to prevent the liquid from splashing and avoid cross-contamination. Close the lid in time after use.

scATAC Barcode Primer use rules/protocols:

Based on the base balancing principle, the scATAC Barcode Primers are used individually or in groups.


- First group: The scATAC Barcode Primer 1 to 4 are used as a base balancing barcode group.
- Second group: The scATAC Barcode Primer 5 to 8 are used as a base balancing barcode group.
- Third group: The scATAC Barcode Primer 9 to 12 are used as a base balancing barcode group.
- Forth group: The scATAC Barcode Primer 13 to 16 are used as a base balancing barcode group.

Totally four groups exist. When all samples have the same data amount requirement, the group of scATAC Barcode Primer varies according to the number of samples. For details, refer to the following table.

Table 30 scATAC Barcode Primer use rules/protocols

Samples/ lane	Method 1	Method 2	Method 3	Method 4
1	1 to 4	5 to 8	9 to 12	13 to 16
2	<ul style="list-style-type: none"> • Sample 1: 1 and 2 • Sample 2: 3 and 4 	<ul style="list-style-type: none"> • Sample 1: 5 and 6 • Sample 2: 7 and 8 	<ul style="list-style-type: none"> • Sample 1: 9 and 10 • Sample 2: 11 and 12 	<ul style="list-style-type: none"> • Sample 1: 13 and 14 • Sample 2: 15 and 16
3	<ul style="list-style-type: none"> • Sample 1: 1 • Sample 2: 2 • Sample 3: 3 and 4 	<ul style="list-style-type: none"> • Sample 1: 5 • Sample 2: 6 • Sample 3: 7 and 8 	<ul style="list-style-type: none"> • Sample 1: 9 • Sample 2: 10 • Sample 3: 11 and 12 	<ul style="list-style-type: none"> • Sample 1: 13 • Sample 2: 14 • Sample 3: 15 and 16
4	<ul style="list-style-type: none"> • Sample 1: 1 • Sample 2: 2 • Sample 3: 3 • Sample 4: 4 	<ul style="list-style-type: none"> • Sample 1: 5 • Sample 2: 6 • Sample 3: 7 • Sample 4: 8 	<ul style="list-style-type: none"> • Sample 1: 9 • Sample 2: 10 • Sample 3: 11 • Sample 4: 12 	<ul style="list-style-type: none"> • Sample 1: 13 • Sample 2: 14 • Sample 3: 15 • Sample 4: 16

Samples/ lane	Method 1	Method 2	Method 3	Method 4
5	<ul style="list-style-type: none"> Sample 1: 1 Sample 2: 2 Sample 3: 3 Sample 4: 4 Sample 5: select any group from the remaining three groups 	<ul style="list-style-type: none"> Sample 1: 5 Sample 2: 6 Sample 3: 7 Sample 4: 8 Sample 5: select any group from the remaining three groups 	<ul style="list-style-type: none"> Sample 1: 9 Sample 2: 10 Sample 3: 11 Sample 4: 12 Sample 5: select any group from the remaining three groups 	<ul style="list-style-type: none"> Sample 1: 13 Sample 2: 14 Sample 3: 15 Sample 4: 16 Sample 5: select any group from the remaining three groups
6	<ul style="list-style-type: none"> Sample 1: 1 Sample 2: 2 Sample 3: 3 Sample 4: 4 Samples 5 and 6: select any two groups from the remaining three groups 	<ul style="list-style-type: none"> Sample 1: 5 Sample 2: 6 Sample 3: 7 Sample 4: 8 Samples 5 and 6: select any two groups from the remaining three groups 	<ul style="list-style-type: none"> Sample 1: 9 Sample 2: 10 Sample 3: 11 Sample 4: 12 Samples 5 and 6: select any two groups from the remaining three groups 	<ul style="list-style-type: none"> Sample 1: 13 Sample 2: 14 Sample 3: 15 Sample 4: 16 Samples 5 and 6: select any two groups from the remaining three groups
7	<ul style="list-style-type: none"> Sample 1: 1 Sample 2: 2 Sample 3: 3 Sample 4: 4 Samples 5 to 7: select groups by referencing the methods used for 3 samples/lane 	<ul style="list-style-type: none"> Sample 1: 5 Sample 2: 6 Sample 3: 7 Sample 4: 8 Samples 5 to 7: select groups by referencing the methods used for 3 samples/lane 	<ul style="list-style-type: none"> Sample 1: 9 Sample 2: 10 Sample 3: 11 Sample 4: 12 Samples 5 to 7: select groups by referencing the methods used for 3 samples/lane 	<ul style="list-style-type: none"> Sample 1: 13 Sample 2: 14 Sample 3: 15 Sample 4: 16 Samples 5 to 7: select groups according to the methods used for 3 samples/lane
8	Select any two groups from the four groups			
8+x (x = 1 to 8, totally 9 to 16 samples)	Perform the following steps: 1. Classify samples 1 to 8 as a group, and add the scATAC Barcode Primer by referencing the methods used for 8 samples/lane. 2. Classify the remaining samples as a group, and correspondingly add different groups of scATAC Barcode Primer based on a value of x by referencing the methods used for 1 to 8 samples/lane.			

 **Tips** It is required to add a mixture of the scATAC Barcode Primer 1 to 16 to a sample. Specifically, take the same volume of the N types of scATAC Barcode Primer, mix them, and add the mixture to the sample.

Appendix 3 96 sequences information and concentration of scATAC Barcode Primer (1 to 96)

This preparation set provides 16 tubes of scATAC Barcode Primer. To meet requirements of multiple samples, this appendix provides 96 sequences information and concentration of scATAC Barcode Primer (1 to 96). You can make the primer according to specific needs.

Preparing scATAC Barcode Primer (1 to 96)

Before preparing primer, use the Nuclease-free Water to dissolve the primer powder into 100 μM solution, and then prepare the primer according to the following table.

 **Tips** You can prepare the primer by percentage according to requirements.

Table 31 scATAC Barcode primer preparation

Name	Volume (μL)
NF Water	80
scATAC Barcode Primer-F (100 μM , NF Water)	10
scATAC Barcode Primer (any one among 1 to 96) (100 μM , NF Water)	10
Total	100

scATAC Barcode Primer (1 to 96) sequences information

Based on the base balancing principle, the scATAC Barcode Primer are used individually or in groups.

- 4 scATAC Barcode Primer (1 to 16) are used as a base balancing barcode group. 1 to 4 are used as a base balancing barcode group, 5 to 8 are used as a base balancing barcode group, and so on.
- 8 scATAC Barcode Primer (17 to 96) are used as a base balancing barcode group. 17 to 24 are used as a base balancing barcode group, 25 to 32 are used as a base balancing barcode group, and so on.
- Details about the sequence information refer to the following table.

 **Tips** The purification method is HPLC.

Table 32 Primer sequence information

Primer	Sequence information (5' → 3')
scATAC Barcode Primer-F	(5'Phosphorylation)GAACGACATGGCTACGATCCGACTT
scATAC Barcode Primer 1	TGTGAGCCAAGGAGTTGATCGGACCTATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 2	TGTGAGCCAAGGAGTTGGATTCCGTCCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 3	TGTGAGCCAAGGAGTTGCGGCAGTAAGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 4	TGTGAGCCAAGGAGTTGTCAATTAGGTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 5	TGTGAGCCAAGGAGTTGCGGATACGAATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 6	TGTGAGCCAAGGAGTTGGCTCGTTACCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 7	TGTGAGCCAAGGAGTTGTTATACGTTGTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 8	TGTGAGCCAAGGAGTTGAACGCGACGTTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 9	TGTGAGCCAAGGAGTTGGCTAGCAGAATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 10	TGTGAGCCAAGGAGTTGCTATCTTCCTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 11	TGTGAGCCAAGGAGTTGAAGCAAGAGCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 12	TGTGAGCCAAGGAGTTGTGCGTGCTTGTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT

Primer	Sequence information (5' → 3')
scATAC Barcode Primer 13	TGTGAGCCAAGGAGTTGCGGATTGCCGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 14	TGTGAGCCAAGGAGTTGGAATCCTGATTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 15	TGTGAGCCAAGGAGTTGTCTGGAATGATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 16	TGTGAGCCAAGGAGTTGATCCAGCATCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 17	TGTGAGCCAAGGAGTTGCATCACTCACTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 18	TGTGAGCCAAGGAGTTGCAGCTGACTCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 19	TGTGAGCCAAGGAGTTGTTGCGAGACATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 20	TGTGAGCCAAGGAGTTGTTGTACCAATTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 21	TGTGAGCCAAGGAGTTGACCACAATCGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 22	TGTGAGCCAAGGAGTTGGGAAGTCTGTTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 23	TGTGAGCCAAGGAGTTGAGAGTGTGGATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 24	TGTGAGCCAAGGAGTTGGCTTGTGGTGTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 25	TGTGAGCCAAGGAGTTGTTGTCCTCTATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 26	TGTGAGCCAAGGAGTTGATTCGCTAGGTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 27	TGTGAGCCAAGGAGTTGCGATGACTACTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 28	TGTGAGCCAAGGAGTTGACAGCTCAGCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 29	TGTGAGCCAAGGAGTTGTATCTAGGTTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT

Primer	Sequence information (5' → 3')
scATAC Barcode Primer 30	TGTGAGCCAAGGAGTTGGAGATGGCAATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 31	TGTGAGCCAAGGAGTTGCGCAAGATCTTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 32	TGTGAGCCAAGGAGTTGGCCGATAGCGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 33	TGTGAGCCAAGGAGTTGCCATCGTTGCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 34	TGTGAGCCAAGGAGTTGTGAACGATTATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 35	TGTGAGCCAAGGAGTTGTAGAGCGAATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 36	TGTGAGCCAAGGAGTTGATGTGTGAGATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 37	TGTGAGCCAAGGAGTTGATCCTAACAGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 38	TGTGAGCCAAGGAGTTGCGCGTCTGCGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 39	TGTGAGCCAAGGAGTTGGATGATCCTTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 40	TGTGAGCCAAGGAGTTGGCTCAACGCTTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 41	TGTGAGCCAAGGAGTTGATGCATCTAATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 42	TGTGAGCCAAGGAGTTGAGCTCTGGACTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 43	TGTGAGCCAAGGAGTTGCTATCACGTGTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 44	TGTGAGCCAAGGAGTTGGGACTAGTGGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 45	TGTGAGCCAAGGAGTTGGCCAAGTCCATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 46	TGTGAGCCAAGGAGTTGCCTGTCAAGCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT

Primer	Sequence information (5' → 3')
scATAC Barcode Primer 47	TGTGAGCCAAGGAGTTGTAGAGGTCTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 48	TGTGAGCCAAGGAGTTGTATGGCAACTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 49	TGTGAGCCAAGGAGTTGCTGCGTACATTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 50	TGTGAGCCAAGGAGTTGATCTCATTAAATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 51	TGTGAGCCAAGGAGTTGAAGTGGCGCATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 52	TGTGAGCCAAGGAGTTGGGCCTTAATGTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 53	TGTGAGCCAAGGAGTTGTCTGAGGCGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 54	TGTGAGCCAAGGAGTTGCGAGCCGATTTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 55	TGTGAGCCAAGGAGTTGGATAACCGGCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 56	TGTGAGCCAAGGAGTTGTCAATATTCCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 57	TGTGAGCCAAGGAGTTGTCCGTTGAATTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 58	TGTGAGCCAAGGAGTTGCAGTACAGTTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 59	TGTGAGCCAAGGAGTTGATTGAGGTACTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 60	TGTGAGCCAAGGAGTTGATTAGAAGTCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 61	TGTGAGCCAAGGAGTTGCAACGCTTCATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 62	TGTGAGCCAAGGAGTTGGGATCGCACGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 63	TGTGAGCCAAGGAGTTGTGCCTCCGATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT

Primer	Sequence information (5' → 3')
scATAC Barcode Primer 64	TGTGAGCCAAGGAGTTGGCGACATCGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 65	TGTGAGCCAAGGAGTTGCATTCTAAGTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 66	TGTGAGCCAAGGAGTTGCAGGCTTGGATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 67	TGTGAGCCAAGGAGTTGATCATCGTCTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 68	TGTGAGCCAAGGAGTTGGTCTTGTGAGTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 69	TGTGAGCCAAGGAGTTGAGTAGGAACGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 70	TGTGAGCCAAGGAGTTGTCACAACCACTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 71	TGTGAGCCAAGGAGTTGGCAGGCCTTCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 72	TGTGAGCCAAGGAGTTGTGGCAAGCTATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 73	TGTGAGCCAAGGAGTTGGAGCATTGTCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 74	TGTGAGCCAAGGAGTTGTGTGATTAGCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 75	TGTGAGCCAAGGAGTTGCCTATGGACTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 76	TGTGAGCCAAGGAGTTGTAGGCGATAGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 77	TGTGAGCCAAGGAGTTGAGACCACGATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 78	TGTGAGCCAAGGAGTTGGTATTAGCCATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 79	TGTGAGCCAAGGAGTTGCTCTGCACTGTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 80	TGTGAGCCAAGGAGTTGACCAGCCTGATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT

Primer	Sequence information (5' → 3')
scATAC Barcode Primer 81	TGTGAGCCAAGGAGTTGGCGTGAGTATTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 82	TGTGAGCCAAGGAGTTGCGCGGAGCATTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 83	TGTGAGCCAAGGAGTTGCAAGTTCACATTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 84	TGTGAGCCAAGGAGTTGAGCACCTCTCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 85	TGTGAGCCAAGGAGTTGTTACAGTGCATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 86	TGTGAGCCAAGGAGTTGTTGCCTAGGCTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 87	TGTGAGCCAAGGAGTTGGCTATGATGGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 88	TGTGAGCCAAGGAGTTGAATTACCATGTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 89	TGTGAGCCAAGGAGTTGAGACATGGTGTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 90	TGTGAGCCAAGGAGTTGCCAGACATATTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 91	TGTGAGCCAAGGAGTTGACGCTTCCTTTTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 92	TGTGAGCCAAGGAGTTGGACGTCTTGATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 93	TGTGAGCCAAGGAGTTGTACTGAGCGGTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 94	TGTGAGCCAAGGAGTTGTGTACACACCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT
scATAC Barcode Primer 95	TGTGAGCCAAGGAGTTGCTTACGTGAATTGTCTTCCTAAGA CCGCTTGGCCTCCGACTT
scATAC Barcode Primer 96	TGTGAGCCAAGGAGTTGGTGTGGAACCTTGTCTTCCTAAG ACCGCTTGGCCTCCGACTT

Appendix 4 Manufacturer information

Manufacturer	Qingdao MGI Tech Co., Ltd
Address	Building 4, No.2, Hengyunshan Road, Qingdao Area, Pilot Free Trade Zone, Shandong, China
Technical support	Qingdao MGI Tech Co., Ltd
Technical support E-mail	MGI-service@mgi-tech.com

---This page is intentionally left blank.---

Part No.: H-020-000746-00

Leading Life Science Innovation

Address: Building 4, No.2, Hengyunshan Road, Qingdao Area,
Pilot Free Trade Zone, Shandong, China

E-mail: MGI-service@mgi-tech.com

Website: www.mgi-tech.com
