

Figure S3: SIRPα-containing complexes are sensitive to Triton X-100 and other detergents. A) Peripheral PBMC, transmigrated monocytes, and BM-derived macrophages (BMDM) were lysed in cold lysis buffer containing 50 mM Tris-HCl, pH 7.4, 150 mM NaCl, 1% Triton X-100, protease inhibitors. After centrifugation, the supernatants were separated and labeled as soluble fractions (S). The pellets were further lysed in lysis buffer containing 0.5% SDS and designated as insoluble fractions (I). Soluble and insoluble fractions were analyzed by SDS-PAGE and WB by anti-SIRPα.ct and anti-Flotillin-1 antibodies. B) Same experiments were performed as in A except using Brij 58 as the detergent.